



www.microflex.com

MICROFLEX GLOVE GUIDE

MICROFLEX
THE MOST TRUSTED NAME IN GLOVES

ISO 13485

P.O. BOX 32000 • RENO, NV 89533-2000
TEL: (800) 876.6866 • TOLL-FREE FAX: (877) 786.4966





THE MOST TRUSTED NAME IN GLOVES®

Microflex® has an innovative approach to glove design and manufacturing. We design gloves to help professionals like you do their job better, and most importantly to protect you when you need it, no matter what the task.

Microflex® gloves provide the comfort, fit and feel that you want, and the high quality, reliable protection you need from The Most Trusted Name in Gloves®.

We know you have many choices when it comes to hand protection. Microflex® brand products provide you with a unique combination: industry-leading quality combined with gloves tailored specifically for your individual needs.

Try Microflex® gloves today and see what makes our gloves stand out from the competition. The difference is quality. The difference is Microflex®.

THE MICROFLEX DIFFERENCE

For over 20 years, Microflex® Corporation has been making gloves specific for your needs. As a leader in the medical exam glove marketplace, Microflex® understands your needs and designs gloves to meet or exceed your expectations.

- *Experience in developing specialized products*
- *Distinct advantage in evolving our product line to meet the needs of the market*
- *Continually improving existing products*
- *On the forefront of new and improved manufacturing methods*

We're proud of our reputation as The Most Trusted Name in Gloves®, and our mission is to do whatever it takes to do better than the rest, to bring you products that consistently deliver the highest quality, comfort and reliability box after box, case after case.

Our mission is our commitment to you and we will continue to work hard every day to achieve the stringent quality and innovation standards we have set.

MICROFLEX AWARDS

2007 Townie Award - Hygienetown

- UltraSense™ won the 2007 "Townie" award (People's Choice award) from Hygienetown
- UltraSense™ was voted the hygienists' favorite non-latex glove in 2007



American Marketing Association 2007

- ACE Award - Print Advertising Campaign - "Microflex Gives You Choices" Campaign
- Award of Excellence - Comprehensive Marketing Program/Campaign - "UltraSense™ re-launch campaign"
- Certificate of Excellence - Print Advertising Campaign - (four color or more) "UltraSense™ re-launch illustrations series"



American Marketing Association 2006

- Certificate of Excellence - Comprehensive Marketing Program / Campaign "The Dawning of a New Era in Gloves"
- Certificate of Excellence - Print Advertising Campaign (four color or more) "The Dawning of a New Era in Gloves"

EMS Expo 2006 - Top 20 Product

- UltraSense™ was awarded the distinction of one of the Top 20 Products from EMS Expo 2006

2006 Frost & Sullivan

- Microflex was awarded Frost & Sullivan's Customer Service Leadership Award for the US disposable medical gloves market
- Per the report, Microflex was awarded this honor due to the introduction of UltraSense™ and ComfortGrip™ in 2006
- This award was given to Microflex for recognizing the need to meet the evolving demands and dynamics of end-users in the healthcare segment by innovating and launching new, value-add products based on ever-changing market conditions in this market



FEATURES AND BENEFITS

- Washed in microprocessor controlled, two-in-one, washer/dryer laundering system to provide very low surface particulates, ionic extractables and non-volatile residues
- Micro textured surface enhances wet or dry grip
- Silicone-free
- Multi-Class cleanroom compatibility
- Provides superb aesthetics, dexterity, tactile sense, fit, feel and comfort
- Hand specific design minimizes unnecessary stress and strain that can cause musculoskeletal injury over time
- Irradiation indicators are placed on all master poly-bags and cartons – your guarantee of sterility assurance

*Note: All poly bags (except carton liner) are manufactured & processed in an ISO 5 cleanroom

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

CE4 12" LATEX
STERILE HAND SPECIFIC GLOVES

ORDER #: HSCE4-879
SIZES: 6.0 - 9.0

GLOVE SPECIFICATIONS

	CE4 12" Latex
Cuff Thickness (mils):	min. 5.1
Palm Thickness (mils):	min. 5.9
Finger Thickness (mils):	min. 7.9
Tensile Strength	
Before aging (MPa):	min. 21
After aging (MPa):	min. 16
Elongation	
Before aging (%):	min. 750
After aging (%):	min. 600
Stress @ 300%	
Before aging:	<2
After aging:	Not applicable

GLOVE ATTRIBUTES

Material:	High grade natural rubber latex
Type:	Sterile
Sizes:	6, 6.5, 7, 7.5, 8, 8.5, 9
Usage:	For single use only
Packaging:	1 pair per poly-wallet, 1 poly-wallet per, vacuum-sealed, inner polybags; 10 inner polybags per, vacuum-sealed, master polybags; 20 master bags per poly lined carton of 200 pairs*



Microflex Glove Products

Quality Assurance	4
Powder-free Synthetic, Standard Exam	5 - 6
Powder-free Latex, Standard Exam	7 - 8
Powder-free Synthetic, Extended Cuff	9 - 10
Powder-free Latex, Extended Cuff	11 - 12
Lightly Powdered Latex, Standard Exam	12 - 14
Lightly Powdered Nitrile, Standard Exam	14
Medical Grade Vinyl, Standard Exam	15
Controlled Environment: Nitrile and Latex	16 - 18
Chemical Resistance Information	18
Marketing Awards 2006	19



CHEMICAL RESISTANCE INFORMATION
AND CUSTOM CHEMICAL TESTING

Microflex® offers a list of chemical ratings for latex, nitrile, and chloroprene glove materials. Ask your authorized distributor representative for a **Microflex® Chemical Resistance Guide**, or visit www.microflex.com for a printable version.



QUALITY IS THE KEY – ISO 13485 CERTIFIED

At the center of our reputation as The Most Trusted Name In Gloves®, is the Microflex® Quality System – industry-leading processes and professionals dedicated to the safety and satisfaction of our valued customers. In addition, Microflex® is THE ONLY glove distributor in the U.S. that is ISO 13485 registered. This quality assurance sets Microflex® apart from the competition.

- *ISO 13485 requirements reach beyond the normal quality system standards of FDA, QSR, and ISO 9000*
- *Distributes into over 35 countries worldwide*
- *Member of the C-TPAT program (Customs Trade Partnership against Terrorism)*
- *Dedicated Complaint Management services for our distributors and end-users*

Microflex® gloves meet or exceed all current and proposed medical examination glove standards imposed by the ASTM International and the Food and Drug Administration (FDA).

CUSTOMS – TRADE PARTNERSHIP AGAINST TERRORISM (C-TPAT)

Microflex® became one of the first companies in the U.S. to join the Department of Homeland Security in a joint government-business initiative to build cooperative relationships that strengthen overall supply chain and border security. Through this initiative, Microflex® has been asked to ensure the integrity of its security practices and communicate security guidelines to its business partners within the supply chain.

Microflex® has also conducted a comprehensive self-assessment of supply chain security using C-TPAT guidelines developed by Customs and the trade community (available for review at www.customs.ustreas.gov). By participating, Microflex® has the opportunity to play an active role in the war against terrorism and ensure greater security for employees, suppliers and customers.

CE4 12" LATEX

AMBIDEXTROUS GLOVES

ORDER #: CE4-200
SIZES: S - XL

GLOVE SPECIFICATIONS

	CE4 12" Latex
Palm Thickness (mils):	6.3 ± 0.8
Finger Thickness (mils):	7.9 ± 0.8
Tensile Strength	
Before aging (MPa):	min. 21
After aging (MPa):	min. 16
Elongation	
Before aging (%):	min. 750
After aging (%):	min. 600
Stress @ 300%	
Before aging:	<2
After aging:	Not Applicable

GLOVE ATTRIBUTES

Material:	Natural rubber latex
Type:	Non-sterile
Sizes:	S - XL
Glove Exterior:	Textured fingertips
Usage:	For single use only
Packaging:	100 gloves per polybag 10 polybags per carton*

*Note: All poly bags (except carton liner) are manufactured & processed in an ISO 5 cleanroom



FEATURES AND BENEFITS

- *Washed in microprocessor controlled, two-in-one, washer/dryer laundering system to provide very low surface particulates, ionic extractables and non-volatile residues*
- *Micro textured surface enhances wet or dry grip*
- *Silicone-free*
- *Multi-Class cleanroom compatibility*
- *Provides superb aesthetics, dexterity, tactile sense, fit, feel and comfort*
- *Double vacuum sealed packaging*

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

CE4 12" NITRILE

AMBIDEXTROUS GLOVES

ORDER #: CE4-313
SIZES: S - XL

GLOVE SPECIFICATIONS

	CE4 12" Nitrile
Palm Thickness (mils):	5.2 ± 1.2
Finger Thickness (mils):	6.7 ± 1.2
Tensile Strength	
Before aging (MPa):	min. 14
After aging (MPa):	min. 14
Elongation	
Before aging (%):	min. 500
After aging (%):	min. 400

GLOVE ATTRIBUTES

Material:	100% nitrile
Type:	Non-sterile
Sizes:	S - XL
Usage:	For single use only
Packaging:	100 gloves per polybag 10 polybags per carton*

NOT MADE FROM NATURAL RUBBER LATEX

*Note: All poly bags (except carton liner) are manufactured & processed in an ISO 5 cleanroom.
U.S. PATENT NO. RE. 35.616. EUROPEAN PATENT NO. 0456333.



FEATURES AND BENEFITS

- *Washed in microprocessor controlled, two-in-one, washer/dryer laundering system to provide very low surface particulates, ionic extractables and non-volatile residues*
- *100% nitrile to alleviate the potential for adverse reactions associated with latex protein*
- *Does not contain plasticizers*
- *Silicone-free*
- *Multi-Class cleanroom compatibility*
- *Provides superb aesthetics, dexterity, tactile sense, fit, feel and comfort*
- *Offers excellent Static Decay and Surface Resistivity properties*
- *Double vacuum sealed packaging*

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



FEATURES AND BENEFITS

- Choose from 9" or 12" cuff lengths to meet your specific protocol
- Low surface particulates, ionic extractables and non-volatile residues
- Micro-textured fingertip surface enhances wet or dry grip
- Silicone free
- Multi-class cleanroom compatibility
- Provides superb dexterity, tactile sense, fit, feel and comfort

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

CE5 9" & 12" LATEX

AMBIDEXTROUS GLOVES

ORDER #: 9" CE5-429; 12" CE5-512
SIZES: S - XL

GLOVE SPECIFICATIONS

	CE5 9/12" Latex
Palm Thickness (mils):	6.3 ± 0.8
Finger Thickness (mils):	7.9 ± 0.8
Tensile Strength	
Before aging (MPa):	min. 21
After aging (MPa):	min. 16
Elongation	
Before aging (%):	min. 750
After aging (%):	min. 600
Stress @ 300%	
Before aging:	<2
After aging:	Not applicable

GLOVE ATTRIBUTES

Material:	Natural rubber latex
Type:	Non-sterile
Sizes:	S - XL
Glove Exterior:	Textured fingertips
Usage:	For single use only
Packaging:	100 gloves per polybag 10 polybags per carton*

*Note: All poly bags (except carton liner) are manufactured & processed in an ISO 5 cleanroom

CE5 9" & 12" NITRILE

AMBIDEXTROUS GLOVES

ORDER #: 9" CE5-637; 12" CE5-755
SIZES: S - XL

GLOVE SPECIFICATIONS

	CE5 9/12" Nitrile
Palm Thickness (mils):	5.2 ± 1.2
Finger Thickness (mils):	6.7 ± 1.2
Tensile Strength	
Before aging (MPa):	min. 14
After aging (MPa):	min. 14
Elongation	
Before aging (%):	min. 500
After aging (%):	min. 400

GLOVE ATTRIBUTES

Material:	100% nitrile
Type:	Non-sterile
Sizes:	S - XL
Usage:	For single use only
Packaging:	100 gloves per polybag 10 polybags per carton*

*Note: All poly bags (except carton liner) are manufactured & processed in an ISO 5 cleanroom

NOT MADE FROM NATURAL RUBBER LATEX

U.S. PATENT NO. RE. 35.616. EUROPEAN PATENT NO. 0456333.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

ULTRASENSE™

POWDER-FREE NITRILE

ORDER #US-220

GLOVE SPECIFICATIONS

	ASTM*	US
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	2.8
Palm Thickness (mils):	min. 2.0	3.5
Finger Thickness (mils):	min. 2.0	4.7
Tensile Strength		
Before aging (MPa):	min. 14	min. 18
After aging (MPa):	min. 14	min. 16
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	400

GLOVE ATTRIBUTES

Material:	Nitrile, not made from natural rubber latex
Shape:	Ambidextrous
Color:	Custom Blue
Sizes:	XS - XL
Cuff:	Beaded
Glove Interior:	Powder-free
Glove Exterior:	Textured fingers
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight; 10 boxes per case



RECOMMENDED FOR USE WHEN:

- Handling small instruments, where fit and feel is key
- Light tasks require excellent grip on a day-to-day basis
- Performing extended or repetitive duties that may cause hand fatigue

Non-stick surface properties resistant to tapes and adhesives

NOT MADE FROM NATURAL RUBBER LATEX

* ASTM D 6319
Patent No. 6,451,893.
Patent No. 7,176,260.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

FREEFORM® SE

POWDER-FREE NITRILE

ORDER #FFS-700

GLOVE SPECIFICATIONS

	ASTM*	FF SE
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	3.5
Palm Thickness (mils):	min. 2.0	3.9
Finger Thickness (mils):	min. 2.0	5.5
Tensile Strength		
Before aging (MPa):	min. 14	min. 18
After aging (MPa):	min. 14	min. 16
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	400

GLOVE ATTRIBUTES

Material:	Nitrile; not made from natural rubber latex
Shape:	Ambidextrous
Color:	Blue
Sizes:	XS - XL
Cuff:	Beaded
Glove Interior:	Powder-free
Glove Exterior:	Textured fingers
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight; 10 boxes per case



RECOMMENDED FOR USE WHEN:

- Working in wet environments that demand a reliable grip
- Light to mid-duty tasks require a combination of comfort and durability
- Precise contact and softness require enhanced tactile sensitivity

NOT MADE FROM NATURAL RUBBER LATEX

* ASTM D 6319
Patent No. 6,451,893.
Patent No. 7,176,260.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



RECOMMENDED FOR USE WHEN:

- *Strength and protection are of great importance*
- *Day-to-day hazards require the best combination of durability and security*
- *Working in environments that demand consistent grip and reliability*

NOT MADE FROM
NATURAL RUBBER LATEX

*ASTM D 6319
U.S. PATENT NO. RE. 35,616. EUROPEAN PATENT NO. 0456333.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

SUPRENO® SE

POWDER-FREE NITRILE
ORDER #SU-690

GLOVE SPECIFICATIONS

	ASTM*	SU SE
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	4.3
Palm Thickness (mils):	min. 2.0	5.5
Finger Thickness (mils):	min. 2.0	7.9
Tensile Strength		
Before aging (MPa):	min. 14	min. 14
After aging (MPa):	min. 14	min. 14
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	400

GLOVE ATTRIBUTES

Material: Nitrile; not made from natural rubber latex
Shape: Ambidextrous
Color: Blue
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Polymer-coated, powder-free
Glove Exterior: Textured fingers
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight, 10 boxes per case

DERMA CARE®

LIGHTLY POWDERED VINYL
ORDER #DC-400

GLOVE SPECIFICATIONS

	ASTM*	DC
Length (mm/in):	min. 230	230/9.1
Cuff Thickness (mils):	n/a	3.5
Palm Thickness (mils):	min. 3.1	4.3
Finger Thickness (mils):	min. 2.0	3.9
Tensile Strength		
Before aging (MPa):	min. 11	min. 11
After aging (MPa):	min. 11	min. 11
Elongation		
Before aging (%):	min. 300	350
After aging (%):	min. 300	300

GLOVE ATTRIBUTES

Material: Vinyl, not made from natural rubber latex
Shape: Ambidextrous
Color: Clear
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Lightly powdered
Glove Exterior: Smooth
Usage: For single use only
Packaging: 100 gloves per dispenser box, by count; 10 boxes per case

FEATURES AND BENEFITS

- *Latex-free formulation*
- *Thicker for added protection without sacrificing tactile sensitivity*
- *Lightly powdered for easy donning*

NOT MADE FROM
NATURAL RUBBER LATEX

*ASTM D 5250

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



RECOMMENDED FOR USE WHEN:

- *Contact with specific** chemicals is a standard activity each day*
- *Tasks require the strength of a synthetic and desire the comfort of latex*
- *Working in wet environments, where an enhanced grip and tactile sensitivity are a must*

NOT MADE FROM
NATURAL RUBBER LATEX

**Tested for use with chemotherapy drugs. Gloves used for protection against chemotherapy drugs exposure must be selected specifically for the type of chemicals used. Tested for use against Carmustine, Cyclophosphamide, Doxorubicin HCl, 5-Fluorouracil, Cisplatin, Etoposide, Paclitaxel, Thio-Tepa, Mitomycin C, Daunorubicin HCl, Ifosfamide, Melphalan, Idarubicin HCl, Carboplatin, Bleomycin Sulfate, Mitoxantrone, Cytarabine HCl.

*ASTM D 6977

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

NEOPRO®

POWDER-FREE CHLOROPRENE
ORDER #NPG-888

GLOVE SPECIFICATIONS

	ASTM*	NP
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	min. 2.0	4.3
Palm Thickness (mils):	min. 2.0	5.1
Finger Thickness (mils):	min. 2.0	6.7
Tensile Strength		
Before aging (MPa):	min. 14	min. 14
After aging (MPa):	min. 14	min. 14
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	400

GLOVE ATTRIBUTES

Material: Polychloroprene
Shape: Ambidextrous
Color: Green
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Polymer-coated, powder-free
Glove Exterior: Textured fingers
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

DERMA FREE™

POWDER-FREE VINYL
ORDER #DF-850

GLOVE SPECIFICATIONS

	ASTM*	DF
Length (mm/in):	min. 230	230/9.05
Cuff Thickness (mils):	n/a	3.5
Palm Thickness (mils):	min. 3.1	4.3
Finger Thickness (mils):	min. 2.0	3.9
Tensile Strength		
Before aging (MPa):	min. 11	min. 11
After aging (MPa):	min. 11	min. 11
Elongation		
Before aging (%):	min. 300	350
After aging (%):	min. 300	300

GLOVE ATTRIBUTES

Material: Vinyl, not made from natural rubber latex
Shape: Ambidextrous
Color: Clear
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Powder-free
Glove Exterior: Smooth
Usage: For single use only
Packaging: 100 gloves per dispenser box, by count; 10 boxes per case

FEATURES AND BENEFITS

- *Latex-free formulation*
- *Thicker for added protection without sacrificing tactile sensitivity*
- *Powder-free*

NOT MADE FROM
NATURAL RUBBER LATEX

*ASTM D 5250

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.





FEATURES AND BENEFITS

- *Low modulus natural rubber latex allows for reduced hand fatigue over long periods of wear*
- *Lightly powdered with USP absorbable cornstarch for easy donning*
- *Longer than typical latex exam gloves*

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

TECH ONE[®]
LIGHTLY POWDERED LATEX
ORDER #MF-101

GLOVE SPECIFICATIONS

	ASTM*	TO
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.9
Finger Thickness (mils):	min. 3.1	6.7
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm ²	max. 150 µg/g

GLOVE ATTRIBUTES

- Material:** Natural rubber latex
- Shape:** Ambidextrous
- Color:** White
- Sizes:** XS - L
- Cuff:** Beaded
- Glove Interior:** Lightly powdered with USP absorbable cornstarch
- Glove Exterior:** Smooth
- Usage:** For single use only
- Packaging:** 100 gloves per dispenser box, by weight; 10 boxes per case

EVOLUTION ONE[®]
POWDER-FREE LATEX
ORDER #EV-2050

GLOVE SPECIFICATIONS

	ASTM*	EV One
Length (mm/in):	min. 220/8.7	255/10.0
Cuff Thickness (mils):	n/a	4.7
Palm Thickness (mils):	min. 3.1	6.7
Finger Thickness (mils):	min. 3.1	8.7
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm ²	max. 50 µg/g

GLOVE ATTRIBUTES

- Material:** Natural rubber latex
- Shape:** Ambidextrous
- Color:** Natural
- Sizes:** XS - XL
- Cuff:** Beaded
- Glove Interior:** Polymer-coated, powder-free
- Glove Exterior:** Fully textured
- Usage:** For single use only
- Packaging:** 100 gloves per dispenser box, by weight; 10 boxes per case

RECOMMENDED FOR USE WHEN:

- *Maximum comfort is required for lasting wear*
- *Rapid donning and anti-sticking properties are needed*
- *Working with a variety of instruments and equipment*

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



FEATURES AND BENEFITS

- *Thickest at the palm and finger for increased protection where you need it most*
- *Fully textured for a secure grip in wet and dry conditions*
- *FOR INDUSTRIAL USE ONLY*

NOT MADE FROM
NATURAL RUBBER LATEX

CAUTION: These gloves are for non-medical use only. They may **NOT** be worn for barrier protection in medical or healthcare applications. Please select other gloves for these applications. Follow your institution's policies for use. For single use only.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

NITRON ONE[®]
LIGHTLY POWDERED NITRILE
ORDER #NO-123

GLOVE SPECIFICATIONS

	NO
Length (mm):	248/9.8
Cuff Thickness (mils):	4.3
Palm Thickness (mils):	5.5
Finger Thickness (mils):	7.9
Tensile Strength	
Before aging (MPa):	min. 14
After aging (MPa):	min. 14
Elongation	
Before aging (%):	500
After aging (%):	400

GLOVE ATTRIBUTES

- Material:** Nitrile; not made from natural rubber latex
- Shape:** Ambidextrous
- Color:** Blue
- Sizes:** S - XL
- Cuff:** Beaded
- Glove Interior:** Lightly powdered with USP absorbable cornstarch
- Glove Exterior:** Fully textured
- Usage:** For non-medical, single use only
- Packaging:** 100 gloves per dispenser box, by weight; 10 boxes per case

DIAMOND GRIP[™]
POWDER-FREE LATEX
ORDER #MF-300

GLOVE SPECIFICATIONS

	ASTM*	DG
Length (mm/in):	min. 220/8.7	245/9.8
Cuff Thickness (mils):	n/a	5.1
Palm Thickness (mils):	min. 3.1	6.3
Finger Thickness (mils):	min. 3.1	9.1
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm ²	max. 50 µg/g

GLOVE ATTRIBUTES

- Material:** Natural rubber latex
- Shape:** Ambidextrous
- Color:** Natural
- Sizes:** XS - XL
- Cuff:** Beaded
- Glove Interior:** Powder-free
- Glove Exterior:** Textured fingers
- Usage:** For single use only
- Packaging:** 100 gloves per dispenser box, by weight; 10 boxes per case

RECOMMENDED FOR USE WHEN:

- *Strength and durability are a requirement on a regular basis*
- *Working with utensils and instruments in mid-to heavy-duty environments*
- *Quality and reliability are what you desire in a hard-working glove*

*ASTM D 3578.

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



RECOMMENDED FOR USE WHEN:

- Durable comfort is required for lasting wear
- Consistent grip is necessary for handling utensils on a day-to-day basis
- Tactile sensitivity is important when working with a number of different instruments

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

DIAMOND GRIP PLUS™

POWDER-FREE LATEX
ORDER #DGP-350

GLOVE SPECIFICATIONS

	ASTM*	DGP
Length (mm/in):	min. 220/8.7	248/9.8
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.1
Finger Thickness (mils):	min. 3.1	6.3
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	700
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 50 µg/g

GLOVE ATTRIBUTES

Material: Natural rubber latex
Shape: Ambidextrous
Color: Natural
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Powder-free, polymer-coated
Glove Exterior: Fully textured
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

COLOR TOUCH®

LIGHTLY POWDERED LATEX
ORDER #CT-133

GLOVE SPECIFICATIONS

	ASTM*	CT
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.9
Finger Thickness (mils):	min. 3.1	6.7
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 150 µg/g

GLOVE ATTRIBUTES

Material: Natural rubber latex
Shape: Ambidextrous
Color: Pink, Green, Blue
Sizes: S (Pink), M (Green), L (Blue)
Cuff: Beaded
Glove Interior: Lightly powdered with USP absorbable cornstarch
Glove Exterior: Fully textured
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

FEATURES AND BENEFITS

- Premium quality, low modulus latex for greater comfort and reduced hand fatigue
- Lightly powdered with USP absorbable cornstarch for easy donning
- Fully textured for a more secure grip
- Meets or exceeds federal guidelines for maximum powder and maximum water extractable protein content
- Lightly scented with peppermint
- Color-coded sizes

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



RECOMMENDED FOR USE WHEN:

- A secure grip is needed in wet and dry environments
- Working with small instruments requires maximum sensitivity
- Comfort, quality, and value are crucial to your needs

*ASTM D 3578.

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

COMFORTGRIP™

POWDER-FREE LATEX
ORDER #CFG-900

GLOVE SPECIFICATIONS

	ASTM*	CFG
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.1
Finger Thickness (mils):	min. 3.1	5.9
Tensile Strength		
Before aging (MPa):	min. 18	min. 18
After aging (MPa):	min. 14	min. 14
Elongation		
Before aging (%):	min. 650	650
After aging (%):	min. 500	500
Protein Rating:	max. 200 µg/dm²	max. 50 µg/g

GLOVE ATTRIBUTES

Material: Natural rubber latex
Shape: Ambidextrous
Color: Natural
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Powder-free, chlorinated, polymer coated
Glove Exterior: Fully textured
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

POWERGRIP®

LIGHTLY POWDERED LATEX
ORDER #PG-199

GLOVE SPECIFICATIONS

	ASTM*	PG
Length (mm/in):	min. 220/8.7	245/9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.9
Finger Thickness (mils):	min. 3.1	6.7
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 150 µg/g

GLOVE ATTRIBUTES

Material: Natural rubber latex
Type: Non-sterile
Shape: Ambidextrous
Color: White
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Lightly powdered with USP absorbable cornstarch
Glove Exterior: Fully textured
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

FEATURES AND BENEFITS

- Low modulus natural rubber latex allows for reduced hand fatigue over long periods of wear
- Lightly powdered with USP absorbable cornstarch for easy donning
- Longer length than typical latex exam gloves
- Fully textured for a secure grip in wet or dry conditions

*ASTM D 3578.

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.





FEATURES AND BENEFITS

- Interior/exterior polymer system provides superior anti-sticking properties resulting in highly effective donning in wet, damp or dry conditions.
- Exceptional durability to handle tough tasks
- Extended cuff for increased protection over the wrist and forearm
- Fine fingertip sensitivity for finding a pulse and starting an I.V.

MEETS NFPA STANDARDS

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

SYNETRON®

POWDER-FREE LATEX
ORDER #SY-911

GLOVE SPECIFICATIONS

	ASTM*	SYN
Length (mm/in):	min. 220/8.7	290/11.4
Cuff Thickness (mils):	n/a	6.7
Palm Thickness (mils):	min. 3.1	9.1
Finger Thickness (mils):	min. 3.1	11.8
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	700
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 50 µg/g

GLOVE ATTRIBUTES

Material: Natural rubber latex
Shape: Ambidextrous
Color: Natural
Sizes: S - XXL
Cuff: Extended, beaded
Glove Interior: Powder-free
Glove Exterior: Fully textured
Usage: For single use only
Packaging: 50 gloves per dispenser box, by weight; 10 boxes per case

ULTRASENSE™ EC

POWDER-FREE NITRILE
ORDER #USE-880

GLOVE SPECIFICATIONS

	ASTM*	US EC
Length (mm/in):	min. 220/8.7	295/11.4
Cuff Thickness (mils):	n/a	3.1
Palm Thickness (mils):	min. 2.0	3.9
Finger Thickness (mils):	min. 2.0	5.1
Tensile Strength		
Before aging (MPa):	min. 14	min. 18
After aging (MPa):	min. 14	min. 16
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	500

GLOVE ATTRIBUTES

Material: Nitrile, not made from natural rubber latex
Shape: Ambidextrous
Color: Custom blue
Sizes: XS - XXL
Cuff: Extended, beaded
Glove Interior: Powder-free
Glove Exterior: Textured fingers
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

* ASTM D 6319
Patent No. 6,451,893.
Patent No. 7,176,260.



RECOMMENDED FOR USE WHEN:

- Strength and protection are of great importance
- Working in environments that demand consistent grip and reliability
- Precise contact requires enhanced tactile sensitivity

Non-stick surface properties resistant to tapes and adhesives

NOT MADE FROM
NATURAL RUBBER LATEX

MEETS NFPA STANDARDS

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



FEATURES AND BENEFITS

- Longer than typical latex examination gloves for added protection where you need it most
- Lightly powdered with USP absorbable cornstarch for easy donning
- Smooth surface for excellent tactile sensitivity
- Meets federal guidelines for maximum powder and maximum water extractable protein content

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

MICRO ONE®

LIGHTLY POWDERED LATEX
ORDER #MO-150

GLOVE SPECIFICATIONS

	ASTM*	MO
Length (mm/in):	min. 220/8.7	245/9.4
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.5
Finger Thickness (mils):	min. 3.1	6.3
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 150 µg/g

GLOVE ATTRIBUTES

Material: Natural rubber latex
Shape: Ambidextrous
Color: White
Sizes: XS - XL
Cuff: Beaded
Glove Interior: Lightly powdered with USP absorbable cornstarch
Glove Exterior: Smooth
Usage: For single use only
Packaging: 100 gloves per dispenser box, by weight; 10 boxes per case

FREEFORM® EC

POWDER-FREE NITRILE
ORDER #FFE-775

GLOVE SPECIFICATIONS

	ASTM*	FF EC
Length (mm/in):	min. 220/8.7	290/11.4
Cuff Thickness (mils):	n/a	4.7
Palm Thickness (mils):	min. 2.0	6.3
Finger Thickness (mils):	min. 2.0	7.5
Tensile Strength		
Before aging (MPa):	min. 14	min. 18
After aging (MPa):	min. 14	min. 16
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	500

GLOVE ATTRIBUTES

Material: Nitrile; not made from natural rubber latex
Shape: Ambidextrous
Color: Blue
Sizes: S - XXL
Cuff: Extended, beaded
Glove Interior: Powder-free
Glove Exterior: Textured fingers
Usage: For single use only
Packaging: 50 gloves per dispenser box, by weight; 10 boxes per case



RECOMMENDED FOR USE WHEN:

- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and security
- Working in environments that demand consistent grip and reliability

NOT MADE FROM
NATURAL RUBBER LATEX

MEETS NFPA STANDARDS

* ASTM D 6319
Patent No. 6,451,893.
Patent No. 7,176,260.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.



RECOMMENDED FOR USE WHEN:

- *Strength and protection are of great importance*
- *Day-to-day hazards require the best combination of durability and security*
- *Working in environments that demand consistent grip and reliability*

NOT MADE FROM
NATURAL RUBBER LATEX

MEETS NFPA STANDARDS

*ASTM D 6319
U.S. PATENT NO. RE. 35.616. EUROPEAN PATENT NO. 0456333.

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

SUPRENO® EC

POWDER-FREE NITRILE
ORDER #SEC-375

GLOVE SPECIFICATIONS

	ASTM*	SEC
Length (mm/in):	min. 220/8.7	295/11.6
Cuff Thickness (mils):	n/a	4.5
Palm Thickness (mils):	min. 2.0	5.5
Finger Thickness (mils):	min. 2.0	7.9
Tensile Strength		
Before aging (MPa):	min. 14	min. 14
After aging (MPa):	min. 14	min. 14
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	500

GLOVE ATTRIBUTES

- Material:** Nitrile; not made from natural rubber latex
- Shape:** Ambidextrous
- Color:** Blue
- Sizes:** S - XXL
- Cuff:** Extended, beaded
- Glove Interior:** Polymer-coated, powder-free
- Glove Exterior:** Textured fingers
- Usage:** For single use only
- Packaging:** 50 gloves per dispenser box, by weight; 10 boxes per case

SAFEGRIP®

POWDER-FREE LATEX
ORDER #SG-375

GLOVE SPECIFICATIONS

	ASTM*	SG
Length (mm/in):	min. 220/8.7	300/11.8
Cuff Thickness (mils):	n/a	7.9
Palm Thickness (mils):	min. 3.1	11.4
Finger Thickness (mils):	min. 3.1	16.5
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 50 µg/g

GLOVE ATTRIBUTES

- Material:** Natural rubber latex
- Shape:** Ambidextrous
- Color:** Blue
- Sizes:** S - XL
- Cuff:** Extended, beaded
- Glove Interior:** Powder-free
- Glove Exterior:** Textured fingertips
- Usage:** For single use only
- Packaging:** 50 gloves per dispenser box, by weight; 10 boxes per case



RECOMMENDED FOR USE WHEN:

- *Added protection is required for extreme environments*
- *Extra coverage over the wrist and forearm is needed for enhanced safety and security*
- *Working in situations requiring intense pressure and strain on the hands*

MEETS NFPA STANDARDS

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



RECOMMENDED FOR USE WHEN:

- *Contact with specific** chemicals is a standard activity each day*
- *Tasks require the strength of a synthetic and the comfort of latex*
- *Working in wet environments, where an enhanced grip and tactile sensitivity are a must*

NOT MADE FROM
NATURAL RUBBER LATEX

MEETS NFPA STANDARDS

**Tested for use with chemotherapy drugs. Gloves used for protection against chemotherapy drugs exposure must be selected specifically for the type of chemicals used. Tested for use against Vincristine, Dacarbazine, Amethopterin Hydrate, Carmustine, Cyclophosphamide, Doxorubicin Hydrochloride, 5-Fluorouracil, Cisplatin, Etoposide, Paclitaxel, Thio-Tepa.

*ASTM D 6977

NEOPRO® EC

POWDER-FREE CHLOROPRENE
ORDER #NEC-288

GLOVE SPECIFICATIONS

	ASTM*	NP EC
Length (mm/in):	220/8.7	300/11.8
Cuff Thickness (mils):	min. 2.0	4.7
Palm Thickness (mils):	min. 2.0	6.3
Finger Thickness (mils):	min. 2.0	8.3
Tensile Strength		
Before aging (MPa):	min. 14	min. 14
After aging (MPa):	min. 14	min. 14
Elongation		
Before aging (%):	min. 500	500
After aging (%):	min. 400	500

GLOVE ATTRIBUTES

- Material:** Polychloroprene
- Shape:** Ambidextrous
- Color:** Green
- Sizes:** S - XXL
- Cuff:** Extended, beaded
- Glove Interior:** Powder-free
- Glove Exterior:** Textured fingers
- Usage:** For single use only
- Packaging:** 50 gloves per dispenser box, by weight; 10 boxes per case

CAUTION: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

ULTRA ONE®

POWDER-FREE LATEX
ORDER #UL-315

GLOVE SPECIFICATIONS

	ASTM*	UO
Length (mm/in):	min. 220/8.7	300/11.8
Cuff Thickness (mils):	n/a	7.1
Palm Thickness (mils):	min. 3.1	9.8
Finger Thickness (mils):	min. 3.1	13.8
Tensile Strength		
Before aging (MPa):	min. 18	min. 21
After aging (MPa):	min. 14	min. 18
Elongation		
Before aging (%):	min. 650	750
After aging (%):	min. 500	600
Protein Rating:	max. 200 µg/dm²	max. 50 µg/g

GLOVE ATTRIBUTES

- Material:** Natural rubber latex
- Shape:** Ambidextrous
- Color:** Natural
- Sizes:** S - XL
- Cuff:** Extended, beaded
- Glove Interior:** Powder-free
- Glove Exterior:** Textured fingers
- Usage:** For single use only
- Packaging:** 50 gloves per dispenser box, by weight; 10 boxes per case



RECOMMENDED FOR USE WHEN:

- *Added protection is required for extreme environments*
- *Extra coverage over the wrist and forearm is needed for enhanced safety and security*
- *Strength is required for your toughest tasks*

MEETS NFPA STANDARDS

*ASTM D 3578

CAUTION: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.